

First aid in the workplace

From eye wash to
wound care



1. **First aid kit**
2. **Eye wash**
3. **Resuscitation**
4. **Wound care**



Further information is available on
www.aluderm.de

A W.Söhngen GmbH brand

For more than 100 years, we have stood for quality, innovation and trust in wound care. For decades, we've been known for our orange first aid kits and a wide range of brands for treating minor and major injuries. But to make sure our age isn't immediately apparent, we've worked on our appearance.

From now on, we're combining our experience and quality in one name and one look:

aluderm⁺

Your benefits

- High recognition value thanks to modern, consistent design
- All important information at a glance
- Innovative product solutions, proven quality
- The new **aluderm⁺** shop contains all updated product information

aluderm⁺ dressings

- For first and second degree burns
- Promote wound healing
- Prevent the formation of scabs
- Reduce the risk of scar formation





First Aid Kit

First aid kits from the brand **aluderm+** in a modern, functional design provide the highest possible level of user friendliness. They can be used with one hand, are clear, stackable, divisible and of course glow in the dark.

Your benefits:

- lights up even without power
- lots of space for the first aid contents
- quick to open
- easy to use
- smart wall mount

Choose the basic contents or expanded contents based on the individual risk of injury in your company. If desired, the security seal can protect the DIN part of the contents from being used too early.

aluderm+ first aid kits are available in two sizes with contents of the usual high quality.

3001407	first aid kit DIN 13157
3001408	first aid kit DIN 13169
3001420	first aid kit DIN 13157 + eye wash
3001440	first aid kit DIN 13157 + metal processing
3001421	first aid kit DIN 13157 + dressings



The security seal protects the contents according to DIN. The additional contents can easily be removed.



DIN 13157 + dressings



DIN 13157 + metal processing

aluderm⁺ eye wash

The sterile aluderm⁺ sodium chloride 0.9% eye wash solution is equivalent to natural tears in the human eye and offers fast relief from irritation caused by dust or dirt. The practical design enables it to be opened easily with just one hand.

2010062 eye wash, 200 ml

2010063 eye wash, 500 ml

aluderm⁺ eye wash refill 5 x 20 ml

The sterile aluderm⁺ sodium chloride 0.9% eye wash solution is equivalent to natural tears in the human eye and offers fast relief from irritation caused by dust or dirt. The ampoules are ideal for workplaces with an increased risk of injury.

2010061 eye wash refill 5 x 20 ml



aluderm⁺ neutralising eye wash, 200 ml

The sterile, aluderm⁺ pH-neutral 4.9% phosphate buffer was developed specially for the rapid neutralisation of acids and alkaline substances in the eyes. In an emergency, the user-friendly bottle with ergonomic eyecup enables effective rinsing of the eye for optimal flow and maximum safety. The carefully considered design enables it to be opened quickly with just one hand. After use, re-flushing with the aluderm⁺ eye wash on the way to see the doctor is recommended.

2010064 neutralising eye wash, 200 ml



aluderm+ eye wash station Basic

With the **aluderm+** eye wash station Basic, you have all eye washes together in one place. This modern system solution has eye washes clearly and within reach in a single, hygienic unit, which can be individually tailored to your needs. The transparent hinged lid ensures a quick overview and immediate access in the event of an emergency. The **aluderm+** eye wash station Basic contains the **aluderm+** eye wash 0.9% sodium chloride 500 ml twice, one **aluderm+** neutralising eye wash, pH neutral 4.9% phosphate 200 ml and two **aluderm+** sodium chloride 0.9% eye wash ampoules. Protected against dust, easy to refill and always ready for use – for maximum safety and efficiency in the workplace.

2010065 eye wash station Basic



aluderm+ eye wash station Solo

The **aluderm+** eye wash station Solo is a compact and user-friendly first aid station for quick and effective aid in the workplace. It has one bottle of **aluderm+** eye wash 0.9% sodium chloride 500 ml and one **aluderm+** sodium chloride 0.9% eye wash ampoule and ensures optimal visibility, easy access and maximum hygiene thanks to its dust-proof housing. Thanks to the slim, space-saving design, the station can be mounted almost anywhere. The contents can be adjusted in a flexible manner, so individual requirements can be met without any problems. Easy to refill and immediately ready for use – the **aluderm+** eye wash station Solo is the ideal solution for workplaces where quick first aid is critical.

2010066 eye wash station Solo



aluderm+ PowerBeat

The AEDs by aluderm+ were developed with a focus on mobility, robustness, user-friendliness and reliability. Each model has been independently tested to ensure that it meets the highest standards. They are equipped with the latest technology and functions including voice instructions, easily understandable graphics, quick analysis for the delivery of the shock and cardiopulmonary resuscitation support. All models have an integrated button for children.

- 0102208 aluderm+ PowerBeat X1
- 0102209 aluderm+ PowerBeat X3
- 0102210V aluderm+ PowerBeat X1
- 0102211V aluderm+ PowerBeat X3


aluderm+ LifePad®

The aluderm+ LifePad® is an auxiliary product for first aiders in the event of cardiac arrest and aims to encourage more people to provide first aid. The step-by-step instructions help them to keep a cool head. The coloured LEDs provide visual feedback to ensure the optimal compression pressure. An acoustic signal shows the optimal compression frequency of 100 beats per minute. The flexible material means the aluderm+ LifePad® can be adapted to any body shape and can be used in people from the age of 12.

The contact surface shows the precise positioning of the hands for resuscitation. Can also be used in combination with a defibrillator (the beep for this can be switched off).



- 0102087 aluderm+ LifePad®




aluderm+ AED wall cabinet

Protect your aluderm+ defibrillator reliably while ensuring rapid access in an emergency. The aluderm+ AED wall cabinet is the ideal storage solution for all common automated external defibrillators (AEDs). Made from robust sheet steel with a powder-coated surface, the cabinet offers optimal protection against dust, moisture and mechanical effects. Thanks to its compact structure, the AED cabinet can easily be mounted in corridors, entrance areas or offices. It contributes to better visibility and quick availability of the defibrillator – a key factor in an emergency.

0102212 AED wall cabinet, L

0102213 AED wall cabinet, S



Optionally available with the



aluderm+ LifePad





aluderm+ wound plaster

The **aluderm+** wound pad is made of a soft non-woven fabric with an aluminium-coated surface that prevents adhesion to the wound and thus enables plasters to be changed with almost no pain. In addition to this, all **aluderm+** wound plasters are suitable for first and second degree burns.

Thanks to their breathable structure, the wound plasters offer a pleasant feeling on the skin with optimal protection against external influences such as dirt and pathogens.

Available as classic, sensitive and waterproof wound plasters.



aluderm+ wound spray

The **aluderm+** wound spray is a local wound rinsing solution for the cleaning, rinsing and moistening of wounds and is suitable for the treatment of abrasions, injuries, cuts and burns. **aluderm+** wound spray is non-irritant, prevents biofilm and can be used from the age of 0.

The alcohol-free composition is particularly gentle on the skin and is therefore also suitable for sensitive skin types, children or sensitive parts of the body. Thanks to the practical spray mechanism, the spray can be applied in a contactless and hygienic manner, reducing the risk of further transmission of pathogens.

The mechanical cleaning of the wound removes foreign bodies to prevent skin irritation and infections. **aluderm+** wound spray thus creates an optimal wound environment in terms of moisture, oxygen and pH to effectively support healing.



First aid materials

The first aid materials by **aluderm+** are characterised by their high-quality processing and innovative composition. Developed for professional use while also being suitable for day-to-day use, they offer reliable assistance for smaller injuries and accidents or acute emergencies.

Whether it's at work, when playing sport or when travelling, **aluderm+** supports you wherever help is needed quickly. The products are well thought through, compact and immediately ready for use – ideal for on the go, in the car, in your handbag or in your rucksack. Thanks to the very good quality and the practical user-friendliness, they can be effortlessly integrated into any emergency situation.

Whether you're on the go on excursions, on holiday or on the road, you need quick access, reliability and hygiene. That's exactly what **aluderm+** offers: easily manageable materials that are hygienically packed and optimally tailored to one another. Small injuries can be treated with care without losing any time and with the good feeling of being equipped as well as you can be.

Dressings and compresses

The **aluderm+** compresses and dressings are crowd-pleasers thanks to their top quality and a special material combination of an aluminium-coated surface, an absorbent core and hydrophobic non-woven PP fabric. The aluminium-coated surface ensures that the wound does not adhere to the wound pad. The wound secretion is quickly transported to the absorbent core, where it remains. This ensures gentle care with minimal pain and the best possible wound healing.

Whether you're treating abrasions, cuts, burns or other skin injuries, **aluderm+** compresses and dressings are designed to offer a safe and comfortable solution both in a professional environment and for day-to-day use.



aluderm⁺
experience trust

aluderm⁺ is a brand of

W.Söhngen GmbH
Platter Straße 84
D-65232 Taunusstein
www.aluderm.de

Version 11-2025

© W.Söhngen GmbH • All rights to product images and design owned by W.Söhngen GmbH.
• Any reproduction, including extracts, requires the permission of W.Söhngen GmbH. Designs
may differ from that shown. The terms and conditions of W.Söhngen GmbH apply.
<https://shop.soehngen.com/agb>